

Table S2

Maker	Clone	Distributor	Catlog	AR	Dilution (Ab)	Incubation time (Ab)	Opal
CD8	C8/144B	Dako	M7103	PH6	1:400	ON	620
Granzyme B (GrB)	GrB-7	Monosan	MON7029-1	PH9	1:200	1H	650
CD3	Polyclonal	Dako	A0452	PH6	1:400	ON	520
Ki-67	MIB-1	Dako	M7240	PH6	1:800	1H	690
FoxP3	236A/E7	AbCam	ab20034	PH9	1:200	ON	540
Cytokeratin (CK)	AE1/AE3	Dako	M3515	PH6	1:1600	1H	570